

EMR-TEK Firewave Operating Manual

The Firewave by EMR-TEK

A compact, high-intensity red and near-infrared (NIR) light lamp designed for personal wellness use.

Product Overview

The Firewave is a portable red and near-infrared light lamp featuring optimized LED drivers, a 236W power output, and a compact design. It delivers low-flicker, high-efficiency light exposure for home use.

Here's how each wavelength works together to provide a complete light experience:

■ Red Light (630 nm & 670 nm) → *Surface-Level Illumination*

- Provides a warm, visible glow that contributes to a refreshed and balanced appearance.
- Enhances the visual brightness and atmosphere of your wellness space.

■ Near-Infrared Light (830 nm) → *Deeper Light Interaction*

- Extends beyond the visible spectrum to create a soothing and immersive light experience.
- Complements red light exposure for a balanced, full-spectrum effect that promotes relaxation and balance.

🔥 Why It's a Game-Changer:

The **Firewave** combines high-powered red and near-infrared light with a compact, portable design—offering a versatile and powerful way to enhance daily wellness routines and create a balanced, revitalizing light experience.

How Firewave Works

Light Experience & Wellness Balance

- **Balanced Light Spectrum** – Red and near-infrared wavelengths work together to create a complete and immersive illumination experience.
- **Natural Harmony** – Designed to emulate the warmth and balance of natural sunlight.
- **Relaxing Ambience** – Encourages a sense of calm and focus through consistent light exposure.

Visual Comfort & Ambient Design

- **Even Light Distribution** – Provides smooth, consistent brightness across the illuminated area.
- **Surface Radiance** – Contributes to a refreshed, balanced appearance in any environment.
- **Comfort-Focused Lighting** – Engineered for a visually soothing and enjoyable experience.

Wellness-Centered Performance

- **Compact & Portable** – Easy to position and move for everyday convenience.
- **High-Efficiency Output** – Delivers steady illumination with minimal flicker for visual comfort.
- **Supports Everyday Balance** – Helps create a calm and restorative atmosphere for daily wellness routines.

Lifestyle Integration

- **Versatile Application** – Ideal for home, office, or relaxation spaces.
- **Ambient Enhancement** – Adds warmth and depth to your environment.
- **Mindful Routine Support** – Complements activities like meditation, stretching, or relaxation.

Technical Specifications

- **Lamp Material:** Metal alloy
 - **Power Output:** 236 W
 - **Actual Power Draw:** 101 W
 - **Wavelengths:** 630 nm, 670 nm, 830 nm
 - **LED Quantity:** 12 LEDs + 200 W chip
 - **Beam Angle:** 60°
 - **Irradiance Levels:**
 - 3 inches: > 200 mW/cm²
 - 6 inches: > 142 mW/cm²
 - **Coverage Area:** 32.7 × 21 cm
 - **Color:** Silver
 - **Modular Connectivity:** No
 - **Modes:**
 - **LED Mode:** Small red and near-infrared LEDs
 - **COB Mode:** Large central red and near-infrared source
 - **LED + COB Mode:** Both active
 - **Mounting:** Pulley system included for hanging on doors
 - **Electrical Requirements:** Built-in voltage converter (use the correct detachable power cord for your region)
 - **Compatibility:** Works with EMR-TEK **Light Stand**
-

Important Safety Precautions

Protective Goggles (Mandatory)

- Always wear the blackout protective goggles provided when using the lamp.
- Do not stare directly at the light.
- If goggles are lost, replacement blackout goggles are available from third-party suppliers.
- Note: Blue-light-blocking glasses do **not** provide protection against red or near-infrared light.

General Warnings

- Not for use by children or infants.
- Do not use if you are photosensitive or are taking photosensitizing medications.
- If you are pregnant or nursing, consult a qualified healthcare professional before use.

- Stop use immediately if discomfort, redness, or irritation occurs.
- Do not operate under the influence of alcohol or substances.
- Do not fall asleep while using the lamp.

Placement & Usage Safety

- LEDs and COB components can become hot — do not touch during or immediately after use.
 - Do not place the lamp on a bed or other soft surface.
 - Indoor use only — do not use outdoors, in a sauna, or in humid/wet environments.
 - Operate in temperatures between -20°C and 40°C , with relative humidity between 45–95%.
 - The **Firewave** is not waterproof — keep away from water and moisture.
-

Operating Instructions

1. Attach the power cord to the lamp and plug it into a wall outlet.
 2. Turn on the power switch on the back of the lamp.
 3. Select your desired mode:
 - **LED Mode:** Small LEDs only
 - **COB Mode:** Large central source only
 - **LED + COB Mode:** Both active
 4. Adjust the session timer with the arrow buttons.
 5. Put on blackout goggles before starting.
 6. Note: Some small LEDs emit near-infrared light and may not appear lit — this is normal.
-

Recommended Session Guidelines

Session Duration

- Begin with ~5 minutes in **LED Mode**.
- Gradually increase to **up to 20 minutes** as you become comfortable.
- Once adjusted, **LED + COB Mode** can also be used for **up to 20 minutes**.

Distance

- Maintain a distance of **12–24 inches**.
- For facial illumination: **LED Mode only** at **18–24 inches**.
- Never position closer than **12 inches**.

Frequency

- Up to **5 times per week** as part of a balanced routine.
- Multiple daily sessions are acceptable if focusing on different areas.

Disclaimer

- The **Firewave** is a general wellness product, not a medical device.
- It is **not intended to diagnose, treat, cure, or prevent** any disease or condition.
- EMR-TEK does not provide medical advice or make medical claims.
- Always follow the recommended usage guidelines and safety precautions.
- EMR-TEK assumes no responsibility for improper use or overuse.

 Tag us using your red light **@emrtekinc**
 Join our VIP affiliate team: **emraffiliateteam@gmail.com**