

EMR-TEK Inferno Operating Manual

The Inferno by EMR-TEK

An Almost Full Body Powerful Red Light Lamp for Energy & Recovery

The Inferno Advantage — Next-Level Wellness Lighting

How the Light Works Together for a Complete Experience

- **Red Light (630nm & 670nm) → Surface-Level Light Exposure**
 - Supports a vibrant, healthy-looking appearance
 - Can help maintain smooth and even-looking skin tone
- **Near-Infrared Light (830nm) → Comfort & Warmth**
 - Offers a deep soothing warmth that feels great after activity
 - Commonly used to support an energized and refreshed sensation

Why It's a Game-Changer:

Inferno delivers a **powerful red & NIR light experience** designed for people who want to:

- ✓ Support an active lifestyle
- ✓ Feel refreshed and ready for the day
- ✓ Enhance relaxation and comfort at home

A professional-grade light experience — now in your wellness space.

How Inferno Works

Energizing & Uplifting Light

- Can help you feel refreshed and ready for the day
- A popular addition to morning or pre-activity routines
- Encourages a natural sense of vitality

Relaxation & Whole-Body Comfort

- The warm near-infrared glow creates a soothing experience
- Ideal after long workdays or workouts
- Complements personal recovery and self-care habits

Appearance & Confidence Boost

- Many users enjoy the healthy-looking glow red light helps support
- Can contribute to a smoother and more radiant visual appearance
- Works well alongside your regular skincare routine

Active Lifestyle Support

- Great for people who want to feel at ease in their body
- A favorite among athletes and anyone who moves a lot
- Helps you unwind and feel good after activity

Technical Specifications

- Lamp Material: Metal alloy
- Power Output: 1400W
- Actual Power Draw: 580W $\pm$ 5W
- Wavelengths: 630 nm, 670 nm, 830 nm
- LED Quantity: 132pcs LED + 200W chip
- Beam Angle: 30°
- Irradiance Levels:
 - 6 inches: >200mW/cm²
 - 12 inches: >70mW/cm²
- Coverage Area: Almost full body
- Color: Silver
- Modular Connectivity: Yes
- The lamp provides **red and near-infrared light exposure** in every operating mode
- Includes a **pulley system** for easy hanging — such as behind a door
- Electrical requirements – the lamp has a built-in voltage converter and can be used in any country. Simply use the correct detachable power cord
- This lamp is compatible with:
 - Heavy Duty Stand
 - **Electric Stand** (*when purchased together with the Inferno*)
 - **Light Stand** (*ideal for horizontal positioning; not recommended for vertical use due to weight*)

Benefits of the Inferno

- The Inferno offers **modular connectivity**, allowing two lamps to be linked and controlled together for a larger light experience.
- It can be paired with any **Inferno** or **Firedragon** model using the built-in **IN/OUT** ports and the included patch cords.

Important Safety Precautions — Please Read & Follow Carefully

Protective Blackout Goggles — Mandatory

- **Always wear the provided blackout protective goggles during every session**
- **Especially important when using light near the face or upper body**
- **Do not look directly at the lamp while it is operating**
- **Replacement blackout goggles are widely available online**
- **Note: Blue-light-blocking glasses do not protect from red or near-infrared light**

User Suitability & Sensitivity

- **Intended for adult use only**
- **Not recommended during pregnancy or while breastfeeding unless approved by a professional**
- **Not recommended for individuals who are sensitive to bright light**
- **If you are unsure how your skin or body may respond, speak with a professional before use**
- **Stop use immediately if you experience discomfort or sensitivity**
- **Follow recommended usage times — excessive exposure is not advised**

Placement & Usage Guidelines

- **Light components may become hot — avoid touching during or immediately after use**
- **Do not fall asleep while using the lamp**
- **Indoor use only — avoid humid spaces, saunas, steamy bathrooms, or outdoor environments**
- **Recommended operating environment: -20°C to 40°C, 45%–95% relative humidity**
- **The Inferno is not waterproof — keep away from water to prevent electrical hazards**

Operating Instructions

- Attach the power cord to the lamp and plug it into a grounded wall outlet
 - Turn on the power switch located on the back of the lamp
 - Select your preferred operating mode(s) using the side panel:
 - **LED Mode** — Small red & near-infrared LEDs only
 - **COB Mode** — Large central red & near-infrared globe only
 - **Both LED & COB Mode** — Full combined light experience
 - Adjust session time using the timer arrows
 - Always wear the provided blackout protective goggles during use
 - **Note:** Not all LEDs appear visibly lit — near-infrared light is invisible to the eye and this is completely normal
 - Choose a time of day that feels comfortable and relaxing for you — many users prefer mornings, but anytime works
-

✓ Session Duration

- Start with **5 minutes in LED Mode only**
 - Gradually increase up to **20 minutes** over several weeks as comfort allows
 - Once accustomed, you may use **COB + LED together**, starting again at **5 minutes** and working up to **20 minutes**
-

✓ Recommended Distance & Frequency

- Maintain a distance of **12–24 inches** from the device
- For the face: **LED Mode only, 18–24 inches** away
- You can align the COB with a chosen area if you'd like to focus the light there
 - Do not position closer than **12 inches**
- Up to **5 sessions per week** recommended to keep a balanced routine
 - Optionally, two daily sessions can be used when focusing on **different body areas**
- Minimal clothing provides more comfortable and consistent light exposure
- In the evening, you can use **LED Mode only** for soft, ambient lighting
 - Simply switch **COB OFF**

Disclaimer

The Inferno is a **general wellness device** intended to support personal comfort and lifestyle routines

- Designed for **home use only**
- This is **not** a medical device
- This product is **not intended** to diagnose, treat, cure, or prevent any disease or condition
- EMR-TEK does not provide medical advice
- Individuals with concerns or unique conditions should speak with a medical professional before use
- Always use as directed and follow all safety precautions
- Overuse or improper operation may lead to unwanted results — the user is responsible for safe use

● **Tag us using your red light @emrtekinc**

📧 **Join our VIP affiliate team:**

emraffiliateteam@gmail.com

