



EMR-TEK®

LED Beauty Mask User Manual



@Emr-tek INC



@emrtekinc



emrprotectionservices@gmail.com



www.emr-tek.com

Thank you for choosing your very own LED Beauty Mask. For your benefits, kindly take a few minutes read the manual before using the product.

Product Overview

This LED Beauty Mask is made with soft, skin-friendly liquid silicone material on both the inner and outer layers. Equipped with 54 groups of advanced 4-in-1 LED bulbs, delivering a total of 216 light beams.

It emits with four wavelengths light at: red 630nm, orange red 610nm, green 520nm, and yellow 580nm,

The hollow out design is aimed to reduces the weight of the mask, eliminating the discomfort of sagging and pressure, making it more comfortable to wear.

This mask is IP67-rated for dust and water resistance, allowing it to be rinsed directly under running water for easy and hygienic cleaning.

Product Specifications



LED Therapy Mask



Controller



Type-C Charging Cable



Adjustable Strap

Material	Liquid Silicone
Net Weight	93g
LED Qty	216pcs
Power	4W
Waterproof	IP67
Input	5V 1A
Battery Capacity	2600 mAh
Wavelength	520nm 580nm 610nm 630nm

LED Light Therapy Benefits

This LED Beauty Mask uses specific wavelengths of light to penetrate deep into the skin, stimulating cellular regeneration, increasing blood circulation, and promoting the production of collagen and elastin. The mask helps to:


- Smooth wrinkles and diminish fine lines
- Detoxify and calm the skin
- Improve skin elasticity
- Brighten and even skin tone
- Soothe sensitive skin
- Speed up wound healing and relieve inflammation
- Reduces hyperpigmentation and age spots

Light Therapy Modes:

- Orange Red Light (610nm) + Red Light (630nm): Stimulates collagen and elastin production, boosts skin elasticity, evens skin tone, reduces fine lines and wrinkles, and repairs pigmentation issues.
- Yellow Light (580nm): Calms sensitive skin, inhibits melanin production, brightens the complexion, enhances skincare absorption, and improves redness and capillary dilation.
- Green Light (520nm): In addition to reducing fine lines and wrinkles green light has also been found to be very useful in treating age spots and hyperpigmentation.

How to Use


1. Controller Overview




2. Charging the Controller

- Before your first use, fully charge the controller.
- During charging, the indicator light will blink. Once all lights are solid and no longer flashing, the controller is fully charged.
- After charging, disconnect the controller.
- Do not use the mask while charging.
- Only use the provided USB cable to avoid malfunction.

3、 Wearing the Mask



Ensure the outer surface of the mask is facing you. Thread the strap with Velcro through the left ear hole and secure it.




Thread the other side of the strap through the mask, adjust to fit comfortably, and secure it.




Wear the mask from the top of your head, adjusting the nose bridge for comfort. If you have a ponytail, thread it through the strap hole.

④



Adjust the strap as needed to ensure the mask fits snugly but comfortably.


4. Connecting the Mask to the Controller



Attach the mask to the controller using the provided cable.


5、Powering On/Off

Press and hold the power button for 1.5 seconds to turn on the device. The display screen will light up.



6、Selecting a Mode

Press the "M" button to switch between light therapy modes: red, Green, or yellow light. Once a mode is selected, the timer will start automatically.



Orange Red + Red Light
610nm+630nm



Green Light
520nm



Yellow Light
580nm

7. Adjusting Energy Levels

Short press the power button to cycle through three energy levels. The higher the level, the brighter the light and the stronger the treatment.

Recommendation:

Begin with the lowest level for your first use. If your skin is not sensitive, gradually increase to higher levels. Sensitive skin should use the lowest setting.




8. Setting Treatment Duration

Press and hold the “M” button for 1.5 seconds to set the treatment duration. Choose from 10, 15, 20, or 30 minutes.

We recommend using the mask for 10 minutes daily, not exceeding a total of 30 minutes per day.

You may also combine treatments, such as 10 minutes of red light followed by 10 minutes of green or yellow light.



Steps for Safe Use

1. Sensitivity Test

Before your first treatment, perform a sensitivity test:

- Place the device on your inner arm.
- Turn on the device and run a 10-minute treatment.
- If you experience discomfort, stop immediately.
- After treatment, wait at least 6 hours. If redness or itching occurs, do not use the mask on your face.

2. Starting Your Treatment

- Cleanse your face and remove all makeup. You may apply toner, but avoid light-sensitive or whitening products.
- Secure the mask snugly on your face.
- Connect the mask to the controller.
- Press and hold the power button for 2 seconds to turn on the device, select your mode, and begin your treatment.
- After treatment, follow your usual skincare routine.
- Recommended use: 3-5 times per week, starting with every other day.

Cleaning and Maintenance

- Clean the mask before and after each use with a soft, dry cloth or wet wipe.
- The mask can be rinsed under running water, but ensure the controller is unplugged and disconnected.
- Do not rub or twist the mask forcefully.
- Store the device in a dry environment, away from direct sunlight.

Safety Precautions

- Not suitable for children, pregnant or breastfeeding women, or individuals with severe melasma or skin conditions.
- Avoid combining strong exfoliants or acne treatments with light therapy.
- Stop use immediately if you experience itching, swelling, or irritation. Consult a healthcare provider if necessary.

Troubleshooting

Issue	Possible Cause	Solution
The mask light is not on or one mode is not working	The USB charging cable is loose	reconnect the USB cable
	Non-matching USB cable	Please use the dedicated USB cable that comes with the mask
	Low battery or depleted	Please charge the controller
Other issues		Please consult the distributor
Note: If the above methods do not resolve the issue, please stop using the device and consult your local distributor.		

Warranty and After-Sales Service

For service, repairs, or inquiries, please contact your local distributor. Warranty does not cover damage caused by improper use.

Repair information

- Customer name _____
- Contact way _____
- Reason for repair _____



EMR-TEK®